

FAT REDUCTION

Muscle Building

One machine, Two Most Popular
Technologies, Building a Perfect
Body




The CoolTrim-EMS combined with Cryo 360° effectively targets the abdominal region by contracting and tightening the muscles, thereby facilitating fat reduction and body sculpting. The Cryo 360° technology freezes fat cells, allowing for their natural elimination. This method provides a safe and efficient means to contour both the face and body without the need for surgical intervention or recovery time, making it particularly suitable for addressing stubborn fat deposits.

SPECIFICATIONS

Specification	Details
Vacuum	0–100 kPa
EMS Frequency	1–150 Hz
Power	5000 W
Intensity	7 Tesla
Screen	12.4 inch
Pulse Width	310 μs ±10%
Cooling	Air Cooling
Voltage	AC 220 V ±10%, 10 A, 50 Hz AC 110 V ±10%, 10 A, 60 Hz
Cryo Temperature	-10°C to 5°C (adjustable)
Sizes of Handle	A Handle: 235×65×47 mm; 235×74×34 mm; 255×65×55 mm; 255×78×40 mm
	B Handle: 180×65×50 mm; 190×65×55 mm; 210×65×70 mm; 220×65×80 mm
	C Handle: 125×50×45 mm; 135×50×50 mm; 145×50×70 mm
	D Handle: 83×30×50 mm

PRINCIPLE OF TREATMENT



BEFORE AND AFTER

