



SONOLIFT PRO

High-Intensity Focused Ultrasound (HIFU) Machine

10D

MICRO PULSE TECHNOLOGY

LIFTING AND FIRMLING COLLAGEN BOOSTER DEVICE



MEDIXOR



FOR LIFTING & TIGHTENING



Forehead
Wrinkle




Crow's Feet




Fine
Wrinkles




Marionette
Lines



Lip
Lines



Cheek
Definition



Neck
Folds

One Machine, Multi-Functional


- IMMEDIATE RESULTS
- NO DOWNTIME
- NO SCARS
- NO SIDE EFFECTS



COLLAGEN SYNTHESIS
DECLINES IN THE DERMIS
OVER TIME.



FOCUSED ULTRASOUND
SHRINKS SKIN BY DELIVERING
ENERGY DIRECTLY INTO IT.




STIMULATED BY FOCUSED
ENERGY, SKIN COLLAGEN
REGENERATES.




BY REGENERATING COLLAGEN,
THE SKIN IS FIRMED AND
TIGHTENED, WHICH ELIMI- NATES
WRINKLES.


For Slimming & Contouring



FAT CELLS AND CELLULITE IS FORMED WITHIN THE BODY.




THE SMASFOCUS X ENERGY DESTROYS THE TARGETED FAT CELLS.



As part of the body's natural healing response, cells regenerate and join together to fill the gap.

THE FAT LAYER IS NATURALLY ELIMINATED FROM THE BODY, RESULTING IN SMOOTH AND WELL-CONTOURED BODY LINES.



Thorax

A xilla

Arms

Flanks

Abdomen

Thighs

Multi-line Handle



Reduced Treatment Time

The groundbreaking S-shaped working mode can reduce nearly twice the treatment time.



Various Cartridges Up to 8 Types

Cartridges 20,000 Shots

Four Cartridges Accompanied by the Machine

FACE



1.5mm, 3.0mm, 4.5mm

BODY (OPTIONAL)



6.0mm, 8.0mm, 10.0mm, 13.0mm, 16.0mm

- **12.0-INCH HD SCREEN**

Ultra-thinborderultra-largedisplay
screen ratio



- **AVIATION-GRADE POWER SUPPLY**

Built-in instant response sensor significantly
enhances safety performance.



01 Mobile Response Sensor

Once the piezoelectric ceramics stop, the handle
immediately shuts down to prevent harm.

02 Allegro MicroSystems Stepper Motor

American imported motor with extended service life.

03 Japanese Piezoelectric Ceramics

Ensures continuous and stable energy output.




Before & After



Before

Af te r



Before

Af te r



Before

Af te r



Before

Af te r



Before

Af te r



Before

Af te r



Before

Af te r



Before

Af te r

Specifications

Specification	Details
Power	1200 W
Frequency	4 MHz - 7 MHz
Screen Width	12.0 Inches
Length	1-2.5 mm Adjustable
Dot Pitch	8-25 mm Adjustable
Energy Output	0.1-2 J/cm ²
Shots output	Zigzag Dotted and Linear Options



MEDIXORA 
PURE AESTHETIC



EMAIL US
sales@medixora.com



OUR OFFICE
Dwarka Sector 24, New Delhi 110077



CALL US
+91-93-190-190-45/46/47/48/49



VISIT OUR WEBSITE
www.medixora.com

